

Dog CD Markers

- Wide range of specificities
- Fluorescent conjugates for flow cytometry
- Dual and triple color reagents
- Manufactured to ISO 9001 standards



Dog CD Markers

AbD Serotec offers a wide range of antibodies raised against dog cell surface antigens.

As a leading supplier of veterinary immunology reagents, we provide an excellent selection of antibodies raised directly against dog CD marker antigens as shown in the table below. Dual and triple color reagents are also included for rapid T lymphocyte profiling by flow cytometry.

Specificity	Clone	Format	Applications	Product Code	Additional Formats
CD3	CA17.2A12	Purified	C, F, IP	MCA1774GA	FITC
CD3/CD4	CA17.2A12, YKIX302.9	FITC + RPE	F	DC046	
CD3/CD8	CA17.2A12, YCATE55.9	FITC + RPE	F	DC047	
CD3/CD4/CD8	CA17.2A12, YKIX302.9, YCATE55.9	FITC + RPE + Alexa Fluor® 647	F	TC014	
CD4	YKIX302.9	Purified	C, F	MCA1038G	Alexa Fluor® 647, FITC, Pacific Blue™
CD4	CA13.1E4	S/N	C, F, IP	MCA1998S	
CD4/CD8	YKIX302.9, YCATE55.9	FITC + RPE	F	DC048	
CD5	YKIX322.3	Purified	C, F, IP	MCA1037G	APC, FITC, RPE
CD8	YCATE55.9	Purified	C, F	MCA1039G	Alexa Fluor® 647, Alexa Fluor® 700, FITC, Pacific Blue™, RPE
CD8 Alpha	CA9.JD3	S/N	C, F, IP	MCA1999S	
CD8 Beta	CA15.4G2	Purified	C, F, IP	MCA1775	
CD11/CD18	YKIX490.6.4	Purified	F, IP	MCA1040G	FITC
CD11a	CA11.4D3	Purified	C, F, IP	MCA1776	
CD11b	CA16.3E10	S/N	C, F, IP	MCA1777S	
CD11c	CA11.6A1	S/N	C, F, IP	MCA1778S	
CD11d	CA11.8H2	Purified	C, F, IP	MCA1779	S/N
CD11d	CA16.3D3	S/N	C, F	MCA2039S	
CD18	CA1.4E9	Purified	C, F, IP	MCA1780	
CD21	CA2.1D6	Purified	C, F, IP	MCA1781R	Alexa Fluor® 647, RPE
CD34	1H6	Purified	F, WB	MCA2411	Alexa Fluor® 647, APC, Biotin, FITC, RPE
CD44	YKIX337.8.7	Purified	C, F	MCA1041G	
CD45	YKIX716.13	Purified	F, IP	MCA1042G	Alexa Fluor® 647, APC, Biotin, FITC, RPE
CD45	CA12.10C12	S/N	C, F	MCA2035S	
CD45R	YKIX753.22.2	Purified	C, CT, F, IP	MCA1043G	FITC, RPE
CD45RA	CA4.1D3	S/N	C, F, WB	MCA2036S	
CD90	YKIX337.217	Purified	F, IP	MCA1036G	
CD107b	AC17	Purified	F, IF, IP, WB	MCA2558	FITC
CDw90	DOG1402	Purified	C, F, WB	MCA1472	
CLAW-H	YKIX312.11.1	Purified	F	MCA1845G	
MHC Class II Monomorphic	YKIX334.2	Purified	C, F, IP	MCA1044G	FITC
MHC Class II Monomorphic	CA2.1C12	S/N	C, F	MCA2037S	

Applications: C = Immunohistology – frozen, CT = Cytotoxic Assays, F = Flow cytometry, IF = Immunofluorescence, IP = Immunoprecipitation, WB = Western Blotting

To browse our entire range of dog antibodies, please visit our website: www.abdserotec.com/VET

North & South America
 Tel: +1 800 265 7376
 Fax: +1 919 878 3751
 Email: sales.us@abdserotec.com

Europe
 Tel: 00800 CALL 4ABD
 (00800 2255 4223)
 Fax: 00800 2FAX 2ABD
 (00800 2329 2223)
 Email: sales.eu@abdserotec.com

Worldwide
 Tel: +44 1865 852 700
 Fax: +44 1865 852 739
 Email: sales@abdserotec.com

www.abdserotec.com

